



**TULSA-PRO<sup>®</sup>**

Inside-Out Prostate Ablation

**Transurethral  
Directional Thermal Ultrasound**

**PROFOUND**


# TULSA-PRO<sup>®</sup>

How we  
see & plan




Real-time  
MRI  
Guidance

How we  
ablate



Transurethral  
Directional Thermal  
Ultrasound

How we  
control




Closed-loop  
Thermal  
Feedback

A prostate solution that is....


Customizable,

1. Whole gland ablation
2. Ablation post radiation failure
3. Targeted ablation
4. Targeted ablation of enlarged prostate with and without malignant lesion




Predictable,

TULSA-PRO uses real-time in-bore MRI for both procedure planning and guidance. The physician draws the precise ablation boundaries during procedure planning. The autonomous robot then predictably rotates and controls the ablation as needed during delivery.



Incision-free.

The physician reaches the prostate through the urethra.






# TULSA-PRO®

How we see & plan




Real-time MRI Guidance

How we ablate



Transurethral Directional Thermal Ultrasound

How we control




Closed-loop Thermal Feedback



## TULSA-PRO® SYSTEM

“Unlike other ultrasound systems on the market, you can monitor the ultrasound ablation process in real time and get immediate MRI feedback of the thermal dose and efficacy.”

Dr. Steven S. Raman  
Professor of Radiology, Urology and Surgery  
David Geffen School of Medicine at UCLA  
Los Angeles, CA



## TULSA PROCEDURE

“The TACT study demonstrated that treatment with TULSA-PRO® provides safe and effective prostate tissue ablation, with little impact on men’s functional ability compared to well-established treatment modalities...”

Dr. Christian Pavlovich  
Professor, Urology and Oncology  
Johns Hopkins University School of Medicine  
Baltimore MD



### Patient Selection

Patients who are MRI cleared are eligible to have the TULSA Procedure.

### Patient Prep & Positioning

While the patient is under anesthesia, MRI guidance is used to ensure the ultrasound applicator and endorectal cooling device are positioned correctly inside the patient.

### Planning

Using intraoperative in-bore MRI, the physician defines safety margins, guides robotic device positioning, and contours the targeted prostate tissue for each transducer element to define the ablation area and volume.

### Delivery

TULSA-PRO combines real-time MRI thermometry and closed-loop control, allowing the physician to actively monitor tissue heating throughout the prostate and nearby critical structures. The system measures the ablation effect in real-time and, through continuous sweeping movement, automatically adjusts the energy delivery.

### Confirmation

After the procedure is complete, contrast-enhanced MRI confirms accurate ablation of prostate tissue.

# TULSA-PRO®

## BENEFITS

### Customizable

#### The flexibility to uniquely tailor the ablation to each patient

Every patient has unique needs; TULSA-PRO allows physicians to customize, so each patient's life does not have to change.

#### The flexibility to ablate prostates of different sizes

TULSA-PRO can be used to ablate both large and small prostates.

#### The flexibility to benefit different types of patients

Schedule a variety of patients in one day.

### Predictable

#### Physician defines the boundary lines and volume to be ablated; the robot follows the instructions

Ablation process is automated, precise, and predictably avoids impacting healthy tissue.

#### Actively protects the urethra and rectum during ablation to preserve men's natural functions

Favorable safety profile demonstrated in Phase I and TACT clinical trial data.

#### Single short procedure with possibility to schedule four per day in routine practice

Actual ablation time is 1-2 cc/min; 2h total procedure time.

### Incision-free

#### Radiation-free

High-intensity directional ultrasound is used to ablate tissue.

#### High throughput

Consistently perform 3-4 procedures in a day.

#### Transurethral inside-out ablation

The entire prostate is accessible (including anterior and posterior).

#### No energy directed through the rectal wall

Inherently safer than outside-in ablation therapy.

#### Patient tolerability

Minimal pain and fast recovery.

#### Cost savings

MR suite is significantly less expensive than an operating room.

#### Corporate Office



Profound Medical Inc.  
2400 Skymark Ave.  
Unit #6 Mississauga, ON L4W 5K5  
Canada

T: +1 647 476 1350  
F: +1 647 847 3739

#### Europe & Middle East Office

Profound Medical GmbH  
Kehrwieder 9  
20547 Hamburg,  
Germany

T: +49 40 8080 93 342  
F: +49 40 8080 93 111



[www.profoundmedical.com](http://www.profoundmedical.com)

[info@profoundmedical.com](mailto:info@profoundmedical.com)

[www.tulsaprocedure.com](http://www.tulsaprocedure.com)

© Profound Medical Inc., 2021.  
TULSA-PRO® is approved for commercial sale in the USA, EU, Canada, and select other countries.  
CAUTION: Rx only device. PROFOUND and TULSA-PRO are trademarks of Profound Medical Inc.  
For more information regarding the regulatory status of the device in other jurisdictions, please contact us at [info@profoundmedical.com](mailto:info@profoundmedical.com)  
108548D

# PROFOUND